


PCR with a smile



Life can be simple
PCR too!

PCR is an invaluable tool in molecular biology research, and at the heart of this application is the DNA polymerase. At Vivantis Technologies, we believe that a successful PCR starts with quality Polymerases. You can choose from a premium selection of our polymerases, for standard PCR or Multiplex PCR, to Hot-Start PCR applications. It is our goal to make PCR a simple and easy process for researchers around the world. With Vivantis Technologies, PCR will be a walk in the park.



The Chromo DNA Polymerase series is a blend of polymerase with inert colour tracer dyes for easy visualization of the addition of polymerase to the reaction and serve as tracking dye during PCR.

Taq DNA Polymerase (recombinant)



Lot :
 Expiry Date :
 Concentration : 5u/μl
 Supplied with : 2ml of 10X ViBuffer A
 1ml of 10X ViBuffer S
 1ml of 50mM MgCl₂

Store at - 20°C

info@vivanttechnologies.com

Description:

Taq DNA Polymerase is a thermostable DNA polymerase. It is suitable for applications requiring high temperature synthesis of DNA. Taq DNA Polymerase catalyzes the polymerization of nucleotides into duplex DNA in the 5' to 3' direction with the presence of Mg²⁺ but maintains the 5' to 3' exonuclease activity.

Features:

- Thermostable enzyme of approximately 94kDa from *Thermus aquaticus*.
- Ultra pure recombinant protein.
- Replicates DNA at 74°C and exhibits a half-life 40 minutes at 95°C.
- Generates mostly 3' dA overhang PCR products which are suitable for TA cloning.

Unit Definition :

1u is defined as amount of enzyme that required to catalyze the incorporation of 10nmoles of dNTP into acid-insoluble material in 30 minutes at 74°C.

Reaction Buffer:

10X ViBuffer A (without MgCl₂):

500mM KCl, 100mM Tris-HCl (pH9.1 at 20°C) and 0.1% Triton™X-100. The buffer is optimized for use with 0.1-0.2mM of each dNTP.

10X ViBuffer S:


160mM (NH₄)₂SO₄, 500mM Tris-HCl (pH 9.2 at 22°C), 17.5mM MgCl₂ and 0.1% Triton™X-100. The buffer is optimized for use with 0.35mM of each dNTP.

Storage Buffer:

20mM Tris-HCl (pH 8.0 at 22°C), 100mM KCl, 0.5% Tween™ 20, 0.5% Nonidet-P40, 0.1mM EDTA, 1mM DTT and 50% glycerol.

Quality Control:

All preparations are assayed for contaminating endonuclease, exonuclease, and non-specific DNase activities. Functionally tested in DNA amplification.



Amplification Using Vivantis Taq DNA Polymerase

Lane M1 : VC 1kb DNA Ladder
 Lane 1.5kb : 1.5kb PCR product generated using 0.2mM dNTPs and 2.0u Vivantis Taq DNA Polymerase.
 Lane 5kb and 8kb : 5kb and 8kb PCR products generated using 0.25mM dNTPs, 2.5u Vivantis Taq DNA Polymerase and 3% formamide.
 Lane 10kb-20kb : 10kb,15kb and 20kb PCR products generated using 0.36mM dNTPs, 2.5u Vivantis Taq DNA Polymerase and 3% formamide.
 Lane M2 : VC Lambda/HindIII Marker

0.5% TAE agarose gel, 5V/cm

SUGGESTED INITIAL PCR CONDITIONS FOR VARIOUS PCR PRODUCT SIZES WITH VIVANTIS DNA POLYMERASES (#PL1201 - 06 / #PL2201 - 06 / #PL3201 - 06 / #PL4201 - 06)
 REACTION MIX (FINAL CONCENTRATION) :

Primers : 0.2 - 1μM Template: Plasmid (0.02-0.2ng) Lambda (0.1 - 150ng) Genomic (0.05-5μg)	Product Size	100bp - 5kb	5kb - 8kb	8kb - 20kb
	dNTP Mix	100μM	200μM	300μM
	ViBuffer (1X)	A	A	S
	Ultrapure DMSO or formamide	-	3%	3%
DNA Polymerase	Refer to below Table (A)			

TABLE (A) : RECOMMENDED UNITS FOR SPECIFIC VIVANTIS DNA POLYMERASES PER 50μL REACTION VOLUME :

Product Size	Taq (#PL1201 - 06)	Max Taq (#PL2201 - 06)	At Taq (#PL3201 - 06)	AtMax Taq (#PL4201 - 06)
0.1 - 5.0kb	2.0	2.0	2.0	2.0
5.0 - 8.0kb	2.5	2.0	2.5	2.0
8.0 - 20.0kb	2.5	2.0	2.5	2.0
>20.0kb	-	-	-	2.0

* This protocol is subjected to changes depending on the template DNA

Product Size	100bp - 5kb	5kb - 8kb	8kb - 20kb
Denaturation	94°C, 2 min	94°C, 2 min	94°C, 2 min
Denaturation	94°C, 30 s	94°C, 12 s	94°C, 12 s
Annealing*	50 - 68°C, 30 s	50 - 68°C, 30 s	50 - 68°C, 30 s
Extension / 1kb	72°C, 30 s	72°C, 45 s	68°C, 1 min
Cycles	25 - 35	25 - 35	25 - 35
Final Extension	72°C, 7 min	72°C, 7 min	68°C, 7 min

* Primer dependant

Product Use Limitation

This product is for research purpose an *in vitro* use only

2X ViRed Taq Master Mix New!


Features

- Suitable for all routine DNA amplification applications
- Stable at 4°C for 6 months, allowing immediate reaction setup without the time-consuming thawing of reagents
- Reduces set-up time and buffer-dye mixing
- Minimizes potential contamination by eliminating several pipetting steps
- Easy confirmation of complete mixing
- No additional loading dye needed – direct loading of final products onto gels
- Generates mostly 3'dA overhang PCR products which are suitable for TA cloning

Ordering information:

Catalogue No.	Description	Pack Size
CLMM01	2X ViRed Taq Master Mix	100 applications


Amplification of 1.5kb DNA fragment from pTZ region DNA using Vivantis 2X ViRed Taq Master Mix



1.0 % TBE agarose gel

- Lane M** : VC 1kb DNA Ladder
- Lane 1** : DNA amplification product generated with 1.25U of Taq DNA Polymerase
- Lane 2** : DNA amplification product generated with 2X ViRed Taq Master Mix (stored at -20°C)
- Lane 3** : DNA amplification product generated with 2X ViRed Taq Master Mix (after 20 freeze-thaw cycles)


Efficiency analysis of Vivantis 2X ViRed Taq Master Mix – minimum & maximum base pair sizes of PCR products generated



1.0 % TBE agarose gel


- Lane M** : VC DNA Ladder Mix
- Lane 1** : DNA amplification product at 100bp generated with 2X ViRed Taq Master Mix
- Lane 2** : DNA amplification product at 1.5kb generated with 2X ViRed Taq Master Mix
- Lane 3** : DNA amplification product at 5kb generated with 2X ViRed Taq Master Mix

Amplification of 1.5kb DNA fragment from pTZ region using Vivantis 2X AtTaq Master Mix



- Lane M** : VC 1kb DNA Ladder
- Lane 1** : DNA amplification product generated with 1.25U of AtTaq DNA Polymerase
- Lane 2** : DNA amplification product generated with 2X AtTaq Master Mix (stored at -20°C)
- Lane 3** : DNA amplification product generated with 2X AtTaq Master Mix after 20 freeze-thaw cycles)

Efficiency analysis of Vivantis 2X AtTaq Master Mix - minimum and maximum base pair size of PCR product generated



- Lane M** : VC DNA Ladder Mix
- Lane 1** : DNA amplification 100bp product generated with 2X AtTaq Master Mix
- Lane 2** : DNA amplification 1.5kb product generated with 2X AtTaq Master Mix
- Lane 3** : DNA amplification 5kb product generated with 2X AtTaq Master Mix

2X AtTaq Master Mix (Hot Start) New!

Features

- Saves time and reduces contamination due to reduced number of pipetting steps
- Stable at 4°C for 6 months, allowing immediate reaction setup without the time consuming thawing of reagents
- Suitable for all routine DNA amplification applications
- Amplification with enhanced specificity, sensitivity and yield
- Amplification with reduced artifacts, such as primer-dimer formation and mispriming in multiplex amplification

Ordering information:

Catalogue No.	Description	Pack Size
PLMM02	2X AtTaq Master Mix	100 applications

2X Taq Master Mix

Features

- Convenient: Ready to Use
- TA cloning Compatible: Generates 3'dA overhangs
- Saves time: Reduced number of pipetting steps
- Stable: Freeze-thaw up to 20 cycles
- Reproductive: Decreases contamination & error rate

Ordering information:


Catalogue No.	Description	Pack Size
PLMM01	2X Taq Master Mix	100 applications

Comparison of Efficiency between Vivantis 2X Taq Master Mix with Supplier F after Freeze-thaw cycles




- Lane M** : VC 1kb DNA Ladder
- Lane 1** : DNA amplification product generated with 2X Taq Master Mix (after 10 freeze-thaw cycles)
- Lane 2** : DNA amplification product generated with 2X Taq Master Mix (after 15 freeze-thaw cycles)
- Lane 3** : DNA amplification product generated with 2X Taq Master Mix (after 20 freeze-thaw cycles)

Amplification of 5Kb DNA Fragment from lambda DNA Using Vivantis 2X Taq Master Mix



- Lane M** : VC 1kb DNA Ladder
- Lane 1** : DNA amplification product generate with 1.25u of Taq DNA Polymerase
- Lane 2** : DNA amplification product generated with 2X Taq Master Mix (stored at -20°C)
- Lane 3** : DNA amplification product generated with 2X Taq Master Mix (after 20 freeze-thaw cycles)

Vivantis
2X ViRed Taq Master Mix
2X AtTaq Master Mix/
2X Taq Master Mix



2X ViRed Taq Master Mix

Product No : CLMM01
Quantity : 100 reactions



Lot :
Expiry Date :
Supplied with : 4 x 625µl 2X ViRed Taq Master Mix*
3ml of Nuclease-free Water
1ml of 50mM MgCl₂

Store at -20°C
*2X ViRed Taq Master Mix consists of Taq DNA Polymerase, Vibuffer A, dNTPs, MgCl₂, inert red dye and stabilizers.



info@vivantechnologies.com

Description :

2X ViRed Taq Master Mix is an optimized ready-to-use 2X concentrated DNA amplification mixture premixed with red color tracking dye. The ViRed Taq Master Mix contains Taq DNA Polymerase, reaction buffer, dNTPs, MgCl₂, inert red dye and stabilizers needed for routine DNA amplification to obtain a wide range of PCR and DNA products up to 8kb. An inert red dye and stabilizers allows direct loading of final products onto gels for electrophoresis. The red color dye migrates at approximately 400bp on 1% agarose in 1X TBE Buffer.

Features:


- Suitable for all routine DNA amplification applications
- Reduces set-up time and buffer-dye mixing
- Minimizes potential contamination by eliminating several pipetting steps
- Easy confirmation of complete mixing
- No additional loading dye needed – direct loading of final products onto gels

Storage and Stability:

- Stable at -20°C for 18 months or at 4°C for 6 months if properly stored
- Stable for 20 freeze-thaw cycles. To avoid frequent freeze-thaw, keeping small aliquots at -20°C is recommended
- For daily use, keeping aliquots at 4°C is recommended

Quality Control:

All preparations are assayed for contaminating endonuclease, exonuclease, and non-specific DNase activities. Functionally tested in DNA amplification.



Amplification of 1.5kb DNA fragment from pTZ DNA region using 2X ViRed Taq Master Mix in a 50µl reaction mixture (1.0% TBE agarose gel).

Lane M : VC 1kb DNA Ladder
Lane 1 : DNA amplification product generated with 1.25u of Taq DNA Polymerase
Lane 2 : DNA amplification product generated with 2X ViRed Taq Master Mix (store at -20°C)
Lane 3 : DNA amplification product generated with 2X ViRed Taq Master Mix (after 20 freeze-thaw cycles)
Lane M : VC 1kb DNA Ladder
Lane 4 : DNA amplification product generated with 1.25u of Taq DNA Polymerase
Lane 5 : DNA amplification product generated with 2X Taq Master Mix (store at -20°C)
Lane 6 : DNA amplification product generated with 2X Taq Master Mix (after 20 freeze-thaw cycles)

RECOMMENDED PROTOCOL FOR 2X ViRed Taq Master Mix:
Gently mix all solutions after thawing. Spin down briefly and keep on ice. Add the following components in a 0.2ml thin walled PCR tube on ice.
For 50µl reaction volume:

Reagent:	Volume	Final Concentration
2X ViRed Taq Master Mix	25µl	*1X
MgCl ₂ (50mM)	Refer to Table (A)	**For more than 1.5mM MgCl ₂
Primers (Fwd / Rev)	Variable	0.1 - 1 µM each
DNA Template	Variable	0.02 - 5µg
Water, nuclease-free	Adjust final volume to 50µl	

**2X ViRed Taq Master Mix contains a fixed final MgCl₂ concentration of 1.5mM. However, higher concentration may be achieved by adding additional MgCl₂. Please refer to Table (A) if higher MgCl₂ concentration is preferred.
Note : Smaller reaction volume may be achieved provided that the same final concentration of each reaction component is maintained.

CYCLING CONDITIONS (100bp-5kb)	
Denaturation	94°C for 2 minutes
Denaturation	94°C for 30 seconds
Annealing	50 - 68°C for 30 seconds
Extension / 1kb	72°C for 30 seconds
Final Extension	72°C for 7 minutes

} 25 - 35 cycles
This protocol may change depending on the template DNA and primers used.

Table (A) : For more than 1.5mM final MgCl₂ concentration

Volume of MgCl ₂ (50mM) stock to add into 50µl reaction mixture (µl)	Final MgCl ₂ concentration (mM)
0.5	2.0
1.0	2.5
1.5	3.0
2.0	3.5
2.5	4.0